

SEQUENCE PARAMETERS

AND

OPTIONS

MODULE FIVE

IMAGE QUALITY

- Image Quality in MRI is a measure of the diagnostic accuracy and appearance of an image.
- It is defined by the contrast of the images, the ability to spatially resolve detail and the signal-to-noise ratio.

CONTRAST DISCRIMINATION

- Contrast is the difference in relative brightness between pixel values and is the result of the signal received from each voxel after scanning.
- The perception of contrast depends on the number of pixel intensities represented by each gray scale.

CONTRAST DISCRIMINATION

- Fourier transform separates the encoded signal into its individual frequency components.
- These signals are translated into a gray scale ranging between 256 shades of gray.


CONTRAST DISCRIMINATION

- Pixels with the most intense signal are assigned the highest value and are the brightest.
- Pixels with the lowest signal intensity are assigned the lowest value and are represented as the darkest.


FACTORS AFFECT CONTRAST

- TR, TE, ETE, ETL, TI
- Pulse Sequence
- Flip Angle
- Relaxation rates
- Hydrogen density
- Flow
- Contrast Media


(TR) TIME OF REPETITION

- TR is the time between successive RF pulses applied to the same slice.
 - TR affects T1 contrast
 - short TR
 - long TR
-  enhances T1 contrast
 minimizes T1 contrast

FSE: TR = 4500MS




FSE: TR = 6500MS



(TE) TIME TO ECHO

- TE is the time from the initial RF pulse to the middle of the listening window (echo).
 - TE affects T2 contrast.
 - short TE
 - long TE
- minimizes T2 contrast
enhances T2 contrast

SE: TE COMPARISON



I

SE


SE

3000 30
0 90
5256 S20
6.0 1A
12:48



01 0003:0020
02 0003:0041

WL
+856
WW
0545



SELECT DNO
0003:0020 -
0003:0041




(TI) TIME OF INVERSION

- TI is the time from the inverting pulse to the start of the 90° or partial RF pulse.
- TI affects T1 contrast.
- short TI*  minimizes T1 contrast
- long TI*  enhances T1 contrast

FSE: FSE - STIR



FSE: FSE - STIR



STIR / FATSAT



PULSE SEQUENCES

- Pulse sequences are a series of RF pulses and gradient applications in an event.
- Pulse sequences are chosen to enhance or accentuate specific tissue contrast.
 - Spin Echo or FSE T1, T2, PD
 - Gradient Echo or FFE T1, T2*, PD
 - IR, STIR T1, Fat Suppression
 - FLAIR Fluid Suppression

FLIP ANGLES

- Flip angles partially flip net longitudinal magnetization into the transverse plane.
- Lowering the flip angles flips less of net longitudinal magnetization into the transverse plane.
- Flip angles vary based on the type of contrast expected.

CONTRAST PARAMETERS

Flip Angle

short

medium

long

Contrast

T2-wt.-like or T2*

PD-wt.

T1-wt.

RELAXATION TIMES

- T1 relaxation times are dependent on the field strength of the magnet.
 - As field strength increases, T1 lengthens.
- T2 relaxation times do not change significantly with magnetic field strength but do change based on chemical and molecular bonding.

RELAXATION TIMES

- Hydrogen Density is the number of hydrogen atoms present in the tissue.
- Hydrogen Density is one of the major determinants of tissue contrast.
- Flow refers to hydrogen in motion during the acquisition. Imaging techniques cause different appearances of flowing spins.
- Contrast Media shortens the T1 & T2 relaxation times of tissues causing a change in contrast.

SPATIAL RESOLUTION

- Is the ability to distinguish or represent small objects defining the sharpness of an image.
- Spatial resolution is controlled by parameters that affect the amount of tissue or signal represented by the pixel or voxel.
- Images have three dimensions:
 - length and width pixel size or resolution
 - depth slice thickness

SPATIAL RESOLUTION FACTORS

- Factors that affect Spatial Resolution
 - voxel size
 - pixel size
 - slice thickness
 - matrix size
 - field of view (FOV)

SPATIAL RESOLUTION

- Anything that makes a geometric change in the size of the image will affect spatial resolution.

SPATIAL RESOLUTION FACTORS

- Voxel Size
 - Voxels are three-dimensional volume elements.
 - pixel resolution x slice thickness = voxel size.
- Pixel Size
 - Pixels are two-dimensional picture elements.
 - pixels have length and width.
- Decreasing the voxel or pixel size increases resolution

SLICE THICKNESS

- Slice thickness
 - As the slice thickness changes so does the amount of information collected within the voxel.
 - The smaller the amount collected the more accurately the pixel is able to represent the information contained in the image.

FSE: TE/SLICE THICKNESS



MATRIX SIZE

- Matrix size
 - Matrix size defines the length and width of the imaged area and is comprised of pixels.
 - Matrix is selected by the number of phase encoding and frequency encoding steps used to create the image.
 - The more steps and the smaller each step the better the spatial resolution.

FIELD OF VIEW (FOV)

- Field of View (FOV)
 - Field of view is defined as the area of interest demonstrated on the image.
 - Field of view is comprised of a matrix with a particular number of pixels and each pixel a specific slice thickness - voxels
 - The smaller the FOV the better the spatial resolution.

$$\text{FOV} = \text{Matrix} \times \text{Resolution (pixel)}$$

SIGNAL-TO-NOISE RATIO

- Signal-to-noise ratio is used to describe the relative contributions to a detected signal made by the true signal and by random superimposed signals or (noise).
- As the signal strength increases the signal-to-noise ratio increases and better image quality results.

SIGNAL-TO-NOISE RATIO

- Noise superimposed on the signal causes the pixel values to oscillate about a mean value.
- Noise blurs the edges of the tissue interfaces reducing edge acuity.
- Noise is created either by the tissues (non-hydrogen substances) or by the imaging system.

SIGNAL-TO-NOISE RATIO FACTORS

- These factors affect signal-to-noise-ratio
 - proton density
 - field homogeneity
 - voxel volume
 - slice thickness, gap
 - TR,TE,TI
 - NEX
 - Field of View

TIME OF REPETITION (TR)

- TR
 - Lengthening the TR allows more time for the tissues to return signal after being disturbed by the RF pulse.
 - Signal-to-noise ratio is increased with increase in TR.

TIME TO ECHO (TE)

- TE
 - The longer the TE the more T2 relaxation has occurred and the more decay of the transverse magnetization.
 - Increasing TE will decrease the S/N ratio.

TIME TO ECHO (TE)

- TE
 - The longer the TE the more T2 relaxation has occurred and the more decay of the transverse magnetization.
 - Increasing TE will decrease the S/N ratio.

NUMBER OF EXCITATIONS (NEX)

- NEX
 - The number of times the information is collected and averaged to create the image, the higher the signal-to-noise ratio.
 - This is also referred to as the number of acquisitions or signal averages.
 - S/N ratio will be improved by the square root of two, when doubling the NEX which is a 41% increase in signal-to-noise.

MATRIX SIZE

- Matrix Size
 - Most matrix size adjustments require that either the pixel size or resolution also change to maintain, decrease, or increase the FOV.
 - Overall, there may be an increased S/N ratio by the sheer fact that the more number of lines within the matrix requires additional phase encoding steps which will increase S/N ratio.

FIELD OF VIEW (FOV)

- Field of View (FOV)
 - The relationship between the matrix and the FOV has a definite impact on the signal-to-noise ratio by determining the size of the pixels represented in the matrix.
 - Usually when the FOV is decreased while maintaining the matrix size, the pixels decrease in size decreasing S/N.

RECEIVER BANDWIDTH

- The receiver bandwidth is the range of frequencies received during the acquisition.
- In general, the narrower the bandwidth, the higher the signal-to-noise ratio.

RECEIVER BANDWIDTH

- The range can be varied in several ways:
 - sampling the same number of data points over a shorter/longer listening window.
 - changing the number of samples or data points while maintaining the listening window.
 - increasing or decreasing the amplitude of the readout gradient.