

Dual-Energy X-ray Absorptiometry (DXA) Capabilities.

LECTURE 9

Learning Objectives

In this presentation we will discuss the capabilities of a Dual-Energy X-ray Absorptiometry (**DXA**) Scanner. This will include the:

- Basic areas that we scan:
 - Lumbar Spine and Hips
- Forearm.
- Total Body Composition Scan.
- Vertebral Fracture Assessment (VFA) AP total Spine.
- Orthopedic Analysis
- Atypical Femur Fracture Scan.
- Lateral lumbar Spine Scan with Bone Mineral Density (BMD) Results.

Dual-Energy X-ray Absorptiometry (DXA)

Is the gold standard in detecting bone density.

It is the gold standard because it covers very large diverse populations and has been around for a long time.

The areas that are generally scanned in a normal clinical setting are the hips and lumbar spine.

From the results of the hip and spine a bone density diagnosis is determined.

AP Spine Scan Positioning

Lie the patient on their back with their head at the right end of the table.

Make sure the patient is straight and lying flat.

Raise the patient's legs and place the positioning block under the legs.

Start the scan in the middle of L5. This is a couple inches below the navel at the top of the iliac crest.


Scan Parameters

Patient Name:	Adria, Sophia A
Patient ID:	9000
Scan Type:	x Lumbar Spine
Scan ID:	A11140201

Scan Length	8.1	in
Scan Width	4.5	in
Line Spacing	0.0397	in
Point Resolution	0.0355	in

2.55 X 0.08 Coll.
140/100 kVp
2.5 mA avg.
15 seconds
60 Hz

**** XRAY DEVICE READY ****
You may begin scanning



Please position the patient for a Lumbar Spine scan

Start Scan Cancel



AP Spine scan

In this scan we want to see all L1-L4.

Start the scan around the middle of L5.

End the scan at the middle of T-12.

Make sure you have the transverse processes of L1-L4.

You should see the T-12 ribs.

Most of the time you should be able to see a small part of the iliac crest.

Referring Physician:




Image not for diagnostic use
k = 1.138, 40 = 48.0
116 x 149

Total

Scan Information:

Scan Date: November 12, 2005 ID: A11120501
 Scan Type: x Lumbar Spine
 Analysis: November 12, 2005 09:48 Version 12.4:3
 Lumbar Spine
 Operator:
 Model: Discovery C (S/N 81202)
 Comment:



DXA Results Summary:

Region	Area (cm ²)	BMC (g)	BMD (g/cm ³)	T-score	PR (%)	Z-score	AM (%)
L1	14.41	14.44	1.002	0.7	108	2.0	129
L2	15.27	16.33	1.069	0.4	104	1.8	123
L3	16.99	19.69	1.159	0.7	107	2.2	127
L4	18.74	21.27	1.135	0.2	102	1.8	121
Total	65.41	71.72	1.096	0.4	105	1.9	124

Total BMD CV 1.0%, ACF = 1.000, BCF = 1.000, TH = 3.855
 WHO Classification: Normal

Ap Spine results

Left side GE Lunar/Right side is Hologic Reports



Region	1		2	
	BMD (g/cm ³)	Young-Adult (%)	T-Score	Age-Matched (%)
L1	0.913	81	-1.8	98
L2	1.023	85	-1.5	102
L3	1.061	88	-1.2	105
L4	1.399	117	1.7	139
L1-L2	0.968	84	-1.5	101
L1-L3	1.003	86	-1.4	103
L1-L4	1.110	94	-0.6	113
L2-L3	1.044	87	-1.3	104
L2-L4	1.169	97	-0.3	116
L3-L4	1.231	103	0.3	122


Scan Parameters

Patient Name: Adria, Sophia A
Patient ID: 9000
Scan Type: x Right Hip
Scan ID: A11140203

Scan Length: in
Scan Width: 4.5 in
Line Spacing: 0.0397 in
Point Resolution: 0.0355 in


2.55 X 0.08 Coll.
140/100 kVp
2.5 mA avg.
11 seconds
60 Hz

**** XRAY DEVICE READY ****
You may begin scanning



Please position the patient for a Right Hip scan

Start Scan Cancel



Please position the patient for a Left Hip scan

Cancel

Hip Scanning Position

Hip Scanning Instructions


Lie the patient on their back on the table with their head at the right end of the table.

Place the **Foot Positioner** under the patient's legs and align the center of the positioner with the table and the patient's midline.

Rotate the patient's entire approximately 25 degrees inward.

Place the medial part of the foot against the positioner, put the strap around the lateral part of the foot, pull the strap tight and Velcro it to the other side of the positioner.

place the cross hair of the laser 3 inches below the greater trochanter and 1 inch medial to the shaft of the femur and start to scan.



Hip Scan

THIS IS WHAT A HIP SCAN LOOKS LIKE ON A HOLOGIC SCANNER (LEFT) AND A GE LUNAR (RIGHT)

An acceptable hip scan includes the following:

- The femoral shaft is parallel to the edge of the scan image.
- The greater trochanter is centered vertically in the window.
- The entire femoral head is visible.
- A 25° internal rotation of the hip showing minimal or no lesser trochanter on the scan image.

After Scanning the Hip.




Image not for diagnostic use
 k = 1.162, d0 = 61.6
 101 x 91
 NECK: -49 x 15



Scan Information:

Scan Date: 23 May 2013 ID: K0523131R
 Scan Type: a Left Hip
 Analysis: 23 May 2013 17:39 Version 12.6:5
 Left Hip
 Operator:
 Model: QDR 4500C (S/N 47998)
 Comment:

DXA Results Summary:

Region	Area (cm ²)	BMC (g)	BMD (g/cm ²)	T - score	PR (%)	Z - score	AM (%)
Neck	4.61	2.87	0.622	-2.0	73	-0.8	87
Troch	10.15	5.18	0.510	-1.9	73	-1.1	82
Inter	17.34	14.83	0.855	-1.6	78	-1.0	85
Total	32.11	22.88	0.712	-1.9	76	-1.0	85
Ward's	1.11	0.50	0.448	-2.4	61	-0.6	87

Total BMD CV 1.0%, ACF = 1.023, BCF = 1.011, TH = 1.136
 WHO Classification: Osteopenia
 Fracture Risk: Increased



Hip Results

LEFT HOLOGIC/RIGHT GE LUNAR

The Fracture Risk Assessment Tool (FRAX) is a fracture risk calculator that estimates an individual's 10-year probability of incurring a hip or other major osteoporotic fracture.

The FRAX[®] tool has been developed to evaluate fracture risk of patients. It is based on individual patient models that integrate the risks associated with clinical risk factors as well as bone mineral density (BMD) at the femoral neck.

Based on FRAX[®] calculations, patients are at high risk for fracture when they have a 3% or higher hip FRAX and a 20% or higher other major osteoporotic fracture.

This tool now comes with every bone density machine and has the ability to apply FRAX when activated.

FRAX Risk Assessment Tool (FRAX[®])

Frax Risk factors

Age	The model accepts ages between 40 and 90 years. If ages below or above are entered, the programme will compute probabilities at 40 and 90 year, respectively.
Sex	Male or female. Enter as appropriate.
Weight	This should be entered in kg.
Height	This should be entered in cm.
Previous fracture	A previous fracture denotes more accurately a previous fracture in adult life occurring spontaneously, or a fracture arising from trauma which, in a healthy individual, would not have resulted in a fracture. Enter yes or no (see also notes on risk factors).
Parent fractured hip	This enquires for a history of hip fracture in the patient's mother or father. Enter yes or no.
Current smoking	Enter yes or no depending on whether the patient currently smokes tobacco (see also notes on risk factors).
Glucocorticoids	Enter yes if the patient is currently exposed to oral glucocorticoids or has been exposed to oral glucocorticoids for more than 3 months at a dose of prednisolone of 5mg daily or more (or equivalent doses of other glucocorticoids) (see also notes on risk factors).
Rheumatoid arthritis	Enter yes where the patient has a confirmed diagnosis of rheumatoid arthritis. Otherwise enter no (see also notes on risk factors).
Secondary osteoporosis	Enter yes if the patient has a disorder strongly associated with osteoporosis. These include type I (insulin dependent) diabetes, osteogenesis imperfecta in adults, untreated long-standing hyperthyroidism, hypogonadism or premature menopause (<45 years), chronic malnutrition, or malabsorption and chronic liver disease
Alcohol 3 or more units/day	Enter yes if the patient takes 3 or more units of alcohol daily. A unit of alcohol varies slightly in different countries from 8-10g of alcohol. This is equivalent to a standard glass of beer (285ml), a single measure of spirits (30ml), a medium-sized glass of wine (120ml), or 1 measure of an aperitif (60ml) (see also notes on risk factors).
Bone mineral density (BMD)	(BMD) Please select the make of DXA scanning equipment used and then enter the actual femoral neck BMD (in g/cm ²). Alternatively, enter the T-score based on the NHANES III female reference data. In patients without a BMD test, the field should be left blank (see also notes on risk factors) (provided by Oregon Osteoporosis Center).

Calculation Tool

Please answer the questions below to calculate the ten year probability of fracture with BM

Country: **UK** Name/ID: [About the risk factors](#)

Questionnaire:

1. Age (between 40-90 years) or Date of birth
 Age: Y: M: D:

2. Sex Male Female

3. Weight (kg)

4. Height (cm)

5. Previous fracture No Yes

6. Parent fractured hip No Yes

7. Current smoking No Yes

8. Glucocorticoids No Yes

9. Rheumatoid arthritis No Yes

10. Secondary osteoporosis No Yes

11. Alcohol 3 or more units per day No Yes

12. Femoral neck BMD (g/cm²)
 T-Score

BMI 23.9

The ten year probability of fracture (%)

Fracture Type	Probability (%)
Major osteoporotic	19
Hip fracture	4.9

[View NOGG Guidance](#)

Calculation Tool

Please answer the questions below to calculate the ten year probability of fracture with BMD.

Country: **US (Caucasian)** Name/ID: [About the risk factors](#)

Questionnaire:

Age (between 40-90 years) or Date of birth
 Age: Y: M: D:

Sex Male Female

Weight (kg)

Height (cm)

Previous fracture No Yes

Parent fractured hip No Yes

Current smoking No Yes

Glucocorticoids No Yes

Rheumatoid arthritis No Yes

10. Secondary osteoporosis No Yes

11. Alcohol 3 or more units per day No Yes

12. Femoral neck BMD (g/cm²)
 Select DXA

Forearm Scan

Reasons a Forearm scan is performed:

- There is more than one standard deviation difference between the hip and the spine scans.
- The degenerative changes in the spine or hips significantly alters the results therefore, that area needs to be thrown out and another area scanned.
- Hardware has been placed in the spine or hips altering the results therefore, those areas need to be disqualified from interpretation and the forearm added.
- The patient is suspected to have hyperparathyroidism.
- The patient exceeds the table weight capacity.

Forearm Positioning



GE Lunar

Scan Parameters

Patient Name:	Adria, Maria
Patient ID:	
Scan Type:	a L.Forearm
Scan ID:	A03160604


Scan Length cm
Scan Width 10.7 cm
Line Spacing 0.1008 cm
Point Resolution 0.0426 cm

6.10 X 0.05 Coll.
140/100 kVp
2.5 mA avg.
31 seconds
60 Hz

**** XRAY DEVICE READY ****
You may begin scanning

Please position the patient for a L.Forearm scan

Start Scan Cancel



Hologic

Forearm Positioning Instruction.

Make sure the forearm is flat and not rotated!!!

Make sure you use the GE Lunar positioner when using a GE Lunar DXA scanner.


Make sure the GE Lunar positioner is correctly positioned. (Refer to the reference manual before you perform this scan!!!)

These can be done sitting in a chair or lying down.


In any case make sure you do it the same every time. This includes using the same chair every time.




Forearm Image

Hologic




GE Lunar





Forearm Results



Total Body Composition Scan

Body Composition by DXA

[John Shepherd](#), [Bennett Ng](#), [Markus Sommer](#), and [Steven B. Heymsfield](#)

Body composition measurements from DXA have been available since DXA technology was developed 30 years ago but are historically underutilized.

DXA is the preferred method for bone and body composition.

The two X-ray attenuations passing through the body can be used to accurately calculate the mass of two different materials given simple algebra and the physical properties of those materials

DXA can measure regional body composition by subdividing the body using specific well-defined cut lines.

Body Composition by DXA

DXA does expose the patient and operator to ionizing radiation, but the dose is very small to both.

Both Hologic and GE systems have the NHANES data integrated into their software to generate Z-scores for various adiposity and lean mass measures.

DXA uses Fat Mass Index (FMI) FMI has a distinct advantage over BMI for defining obesity status since it is independent of lean mass status.

sex-specific cut points for normal, excess fat, obesity I, obesity II, and obesity III. DXA system also report estimates for visceral adipose tissue (VAT) for either Hologic or GE systems.

Body Composition by DXA

Much of the advancement in the use of DXA for sports and fitness has been application of the technology that has been available for years. Yet, there is a unique place for DXA in evaluating the success of sports, diet, and fitness interventions because of its unique ability to simultaneously measure bone, lean, and fat mass status.

Total body Scan Prep

Restrictions

- No Barium X-ray or Nuclear medicine test within 2 weeks of the study
- No CT (Iodine) or MRI (Gadolinium) contrast within 1 week of the study
- DXA is an x-ray source so we can scan individuals with metal implants. However, large amounts of implants like Rods in the limbs and spine and large plates and screws can affect an accurate outcome.

Day of Scan

Come fasting for six hours

Come not having exercised for 6 hours

Come with having removed ALL body Piercings.

Subject will fill out all appropriate paper work if any.

Subject will be asked to void Bladder.

Subject will Change into scrubs. Subject will need to remove all clothing including socks, for women Bra's will also need to be removed. Underwear can stay on.

Subject will get total body composition scan.

Subject can change back into civilian cloth

Females of reproductive age (14-65) will need a UA (Pregnancy Test) Before the Scan.

Scan parameters

Scan will take approximately 3 minutes.

Patient will be asked to lay on back, with no pillow, with legs and arms straight.

Toes will be pointed up and slightly rotated inward (Hologic). We will have strap to hold feet into position.

If patient moves the scan will have to be started over.

Patient will be asked to not talk during scan.

Scan Parameters

Patient Name: Adria, Sophia A
Patient ID: 9000
Scan Type: a Whole Body
Scan ID: A10180207

Scan Length in
Scan Width in
Line Spacing in
Point Resolution in

4.00 X 0.08 Coll.
140/100 kVp
0.8 mA avg.
165 seconds
60 Hz


**** XRAY DEVICE READY ****

You may begin scanning



Please position the patient for a Whole Body scan

The Regions are predefined as a template overlay and can be adjusted to match the patient's anatomy.



Total Body Scans Analysis Page

Whole Body Fan Beam Analysis

TBAR1058 - NHANES BCA
 F.S. 68.00% -10.00%
 Head assumes 17.0% brain fat
 LBM 73.2% water

Region	Fat (g)	Lean+BMC (g)	Total (g)	%Fat
L Arm	948.5	1662.4	2610.9	36.3
R Arm	717.7	1481.2	2198.9	32.6
Trunk	6689.1	18605.6	25294.7	26.4
L Leg	2764.1	8552.4	11316.5	24.4
R Leg	2899.8	8169.8	11069.5	26.2
Sub Tot	14019.1	38471.4	52490.5	26.7
Head	879.6	2940.3	3819.9	23.0
TOTAL	14898.7	41411.7	56310.4	26.5
Android	1463.8	3181.3	4645.1	31.5
Gynoid	2075.3	6244.4	8319.7	24.9

QDR Workstation SN: 123456
 Version 13.0 :3 05/22/2008 15:26

Results button

Hologic QDR Workstation

Regions A/G Regions Results Sub Regions Subregion Results

Results Toolbox

BMD BCA Rulers

Display Close

Dual Energy

For Help, press F1

Patient Data

Scan Date: 12/21/2009 13:37
 Name: Adria, Ada
 Pat Id: Sex: F
 Birthdate: 07/01/1950 Age: 59
 Height: 67.0 in Weight: 110.0 lb
 Ethnic: White

Auto Whole Body Fan Beam Analysis

TOTAL BMC and BMD CV is < 1.0%

Region	C.F.	Area(cm ²)	BMC(g)	BMD(g/cm ³)
L Arm	1.025	193.89	146.05	0.753
R Arm	1.006	216.24	172.90	0.800
L Ribs	1.000	110.91	86.32	0.778
R Ribs		116.50	80.36	0.690
T Spine		154.00	155.78	1.012
L Spine		60.24	73.25	1.216
Pelvis		194.29	245.12	1.262
L Leg		367.44	430.73	1.172
R Leg		373.43	439.20	1.176
Sub Tot		1786.94	1829.70	1.024
Head		216.24	603.37	2.790
TOTAL		2003.17	2433.08	1.215

Total T-score: 1.3 Total Z-score: 1.5

QDR Workstation SN: 1234
 Version 13.1 :3 12/21/2009 13:39


Help 12/21/2009 13:39

Analysis Results



DXA Results Summary:

Region	BMC (g)	Fat Mass (g)	Lean Mass (g)	Lean+ BMC (g)	Total Mass (g)	% Fat
L Arm	142.88	1723.6	1711.5	1854.4	3578.0	48.2
R Arm	162.16	1559.1	1764.3	1926.4	3485.6	44.7
Trunk	576.13	13157.8	24272.5	24848.7	38006.5	34.6
L Leg	407.72	4285.5	5929.1	6336.8	10622.3	40.3
R Leg	397.16	4239.0	5913.2	6310.4	10549.4	40.2
Subtotal	1686.05	24965.1	39590.6	41276.7	66241.8	37.7
Head	444.22	1033.2	2759.6	3203.9	4237.0	24.4
Total	2130.27	25998.3	42350.3	44480.5	70478.8	36.9

Compartmental Trending



Total Body % Fat



Total Body % Fat Results

Scan Date	Age	% Fat	Change Previous
24.09.2016	29	41.0	-6.6
31.05.2016	29	47.6	

Total Fat Mass Results

Scan Date	Age	Fat Mass (g)	Change Previous
24.09.2016	29	33358	-14852
31.05.2016	29	48210	

Total Lean Mass Results

Scan Date	Age	Lean (g)	Change Previous
24.09.2016	29	45363	-5088
31.05.2016	29	50451	

Total Mass Results

Scan Date	Age	Mass (g)	Change Previous
24.09.2016	29	81351	-19869
31.05.2016	29	101220	

Vertebral Fracture Assessment (VFA)

VFA is an image of the lumbar and thoracic spine acquired on DXA scanners, for the purpose of diagnosing Osteoporotic Vertebral Fractures.

VFA provides the ability to analyze the shape of vertebra using vertebral height measurements and calculations of those heights to compute wedge ratios

VFA is also called IVA by Hologic.

Reasoning for Vertebral Fracture Assessment (VFA)

Two thirds of fractures on X-ray are unappreciated and only 10% necessitated admission to the hospital.

Spinal X-rays may not be routinely performed in osteoporosis evaluation.

Many patients with vertebral fractures are osteopenic by DXA T-score. VFA would allow appropriate diagnosis of osteoporosis in these patients.

Many studies document association between prevalent vertebral fractures and the risk of future vertebral and non-vertebral fractures.

Other Names



Instant Vertebral Assessment

Lateral Vertebral Assessment

Morphometry of Vertebrae

Detection of Fractures

Separate read from DXA

Can eliminate additional spine films

Vertebral Fracture Assessment (VFA)

The thoracolumbar spine is imaged to detect vertebral fractures using DXA.

Lateral images are obtained at time of DXA testing.

Allows integration of BMD with presence of vertebral fracture for diagnosis, fracture risk assessment and patient management.



The screenshot displays the Hologic QDR Workstation software interface for Vertebral Fracture Assessment (VFA). The main window shows a lateral DXA image of the spine. The interface includes a toolbar with buttons for 'DE Scan', 'Multi View', 'Visual Tools', and 'Analysis Tools'. The 'Analysis Tools' section contains a 'L' field with the value '1944', a 'W' field with the value '1156', and buttons for 'Revert', 'Invert', and 'Flip'. There are also zoom in and zoom out icons. The 'Print Report', 'Print Image', and 'Close' buttons are located at the bottom of the toolbar. The patient data section on the right includes fields for Name, Patient ID, Sex, Birth Date, Physician, Height, Weight, Age, Ethnicity, and Scan Date. The vertebral assessment table is also visible, with columns for 'Label', 'Post', 'Mid', 'Ant', 'Wedge', 'Biconcave', and 'Crush' Deformation (Grade). The table is currently empty. The status bar at the bottom indicates 'For Help, press F1' and the date/time '07/17/2013 17:48'.

DE Scan
Multi View
Visual Tools | Analysis Tools
L 1944
W 1156
Revert
Invert
Flip
Print Report
Print Image
Close

Hologic QDR Workstation - [Analyze]
A07171306 a SE Lateral Image Patient Data Deformity Tools

Patient Data
Name: IVA
Patient ID:
Sex: Female
Birth Date: 01/01/1950
Physician:
Height: 65.0 in
Weight: 150.0 lb
Age: 63
Ethnicity: White
Scan Date: 07/17/2013

Vertebral Assessment
Height(mm)
Label Post Mid Ant Wedge Biconcave Crush
Deformity (Grade)

* denotes automatic assessment update
Std Dev 1.0 1.0 1.0 5.0% 5.0% 5.0%

For Help, press F1 Help 07/17/2013 17:48

Lateral Positioning Hologic



C & W Models



SL & A Models
(Rotating Gantry)

Hologic ODR Workstation - [Analyze]

07/17/2013 a SE AP Image Patient Data Deformity Tools

DE Scan
Multi View

Visual Tools Analysis Tools
L 1945

W 1145

Reset
Invert
Flip

Print Report
Print Image
Close



* denotes automatic assessment update

Patient Data

Name: IVA
 Patient ID: 9000
 Sex: Female
 Birth Date: 01/01/1950
 Physician:
 Height: 65.0 in
 Weight: 150.0 lb
 Age: 63
 Ethnicity: White
 Scan Date: 07/17/2013

Vertebral Assessment

Label	Height(mm)			Percent Deformation		
	Post	Mid	Ant	Wedge	Biconcave	Crush
	1.0	1.0	1.0	5.0%	5.0%	5.0%

Std Dev 1.0 1.0 1.0 5.0% 5.0% 5.0%

Help 07/17/2013 17:32

For Help, press F1


Scan Parameters

Patient Name: Adria, Sophia A
 Patient ID: 9000
 Scan Type: f SE AP Image
 Scan ID: A01030305

Scan Length: 16.1 in
 Scan Width: 4.5 in
 Line Spacing: 0.0142 in
 Point Resolution: 0.0355 in

2.40 X 0.02 Coll.
 140/140 kVp
 5.0 mA avg.
 10 seconds
 60 Hz


**** XRAY DEVICE READY ****
 You may begin scanning


Please position the patient for a SE AP Image scan

Start Scan Cancel

LVA Morphometry




Region ¹	Avg. Ht. ²		M/P Ratio ²		A/P Ratio ²	
	(cm)	Z-score	(%)	Z-score	(%)	Z-score
T4	1.70	-1.6	100	1.1	101	0.9
T5	1.97	-0.2	89	-0.5	94	-0.1
T6	1.71	-1.9	99	1.0	95	0.5
T7	1.82	-1.5	80	-1.9	62	-4.5
T8	1.93	-1.1	88	-0.6	59	-4.9
T9	2.04	-0.9	78	-2.3	93	0.0
T10	2.20	-0.8	87	-1.1	92	-0.3
T11	2.37	-0.4	88	-0.7	91	-0.2
T12	2.60	0.0	92	-0.2	97	0.7
L1	2.77	0.1	88	-0.8	99	1.1
L2	2.83	0.0	90	-0.5	93	-0.3
L3	2.80	-0.4	91	-0.5	99	0.2
L4	2.97	0.4	102	0.8	108	1.2

Severe Wedge
 Mild Biconcavity

GE Lunar Report

COMMENTS:




Normal (Grade 0)	<input type="checkbox"/> No Vertebral Deformities Seen		
	Wedge Deformity	Biconcave Deformity	Crush Deformity
Mild (Grade 1)			
Moderate (Grade 2)			
Severe (Grade 3)			

Label	Height (mm)			Percent Deformation		
	Post	Mid	Ant	Wedge	Biconcave	Crush
	Deformity (Grade)					
T4	17.8	18.3	18.2	0.0%	-3.1%	2.5%
T5	19.9	17.9	18.5	6.8%	9.9%	0.0%
T6	19.4	18.1	19.2	1.1%	6.7%	0.0%
T7	21.1	19.6	18.5	12.3%	6.8%	0.0%
T8	20.9	19.4	18.7	10.5%	7.0%	0.0%
T9	21.1	20.2	19.5	7.3%	4.3%	0.0%
T10	22.4	20.2	18.3	18.5%	9.9%	0.0%
T11	24.3	18.8	16.1	33.7%	22.5%	0.0%
T12	26.3	19.4	19.4	26.2%	26.5%	0.0%
L1	28.9	22.5	18.1	37.6%	22.4%	0.0%
L2	27.0	23.7	27.1	0.0%	12.2%	0.5%
L3	27.8	25.4	26.4	5.2%	8.8%	0.0%
L4	27.3	25.1	26.7	2.2%	8.3%	0.0%
Std Dev	1.0	1.0	1.0	5.0%	5.0%	5.0%

Hologic Report

Name: **Physician's Report Writer Rate of Change** Sex: Female Height: 62.0 in
 Patient ID: F219 Ethnicity: White Weight: 140.0 lb
 Age: 54 Date of Birth: 08.07.1946



For Vertebral Deformity Evaluation Only
 Scan Date: June 20, 2001
 Scan ID: A06200108

Scan Type: fSEML Lateral Image




For Vertebral Deformity Evaluation Only
 Scan Date: June 20, 2001
 Scan ID: A06200105


Scan Type: fSEAP Image




Scan Date: June 20, 2001
 Scan ID: A06200106



Scan Type: fLeft Hip



Scan Date: June 20, 2001
 Scan ID: A06200104



Scan Type: fLumbar Spine

Results:

	BMD (g/cm ³)	T-Score	FR (%)	Z-Score	AM (%)	Classification
Left Hip (Neck)	0.734	-1.0	86	0.0	100	Normal
Left Hip (Total)	0.917	-0.2	97	0.5	107	Normal
Spine (Total)	1.023	-0.2	98	0.8	110	Normal

Total BMD CV1%.

Summary:

	Classification
Left Hip BMD (Neck)	Normal
Left Hip BMD (Total)	Normal
Spine BMD (Total)	Normal

A spine fracture indicates 5X risk for subsequent spine fracture and 2X risk for subsequent hip fracture.

World Health Organization criteria for BMD interpretation classify patients as Normal (T-score at or above -1.0), Osteopenic (T-score between -1.0 and -2.5), or Osteoporotic (T-score at or below -2.5).



New Formats



Hologic QDR Workstation - [Analyze]

AD7171306 a SE Lateral Image

DE Scan Multi View

Visual Tools | Analysis Tools

L 1344

W 1156

Revert Invert Flip

Print Report Print Image Close

* denotes automatic assessment update

Patient Data Deformity Tools

Patient Data


Name: N/A
 Patient ID:
 Sex: Female
 Birth Date: 01/01/1950
 Physician:
 Height: 65.0 in
 Weight: 150.0 lb
 Age: 63
 Ethnicity: White
 Scan Date: 07/17/2013


Label	Vertebral Assessment Height(mm)			Percent Deformation		
	Post	Mid	Ant	Wedge	Biconcave	Crush
	1.0	1.0	1.0	5.0%	5.0%	5.0%

Std Dev 1.0 1.0 1.0 5.0% 5.0% 5.0%

Help 07/17/2013 17:48

For Help, press F1






Region ¹	Avg. Ht. ²		M/P Ratio ²		A/P Ratio ²	
	(cm)	Z-score	(%)	Z-score	(%)	Z-score
T4	1.70	-1.6	100	1.1	101	0.9
T5	1.97	-0.2	89	-0.5	94	-0.1
T6	1.71	-1.9	99	1.0	95	0.5
☞ T7	1.82	-1.5	80	-1.9	62	-4.5
☞ T8	1.93	-1.1	88	-0.6	59	-4.9
☞ T9	2.04	-0.9	78	-2.3	93	0.0
T10	2.20	-0.8	87	-1.1	92	-0.3
T11	2.37	-0.4	88	-0.7	91	-0.2
T12	2.60	0.0	92	-0.2	97	0.7
L1	2.77	0.1	88	-0.8	99	1.1
L2	2.83	0.0	90	-0.5	93	-0.3
L3	2.80	-0.4	91	-0.5	99	0.2
L4	2.97	0.4	102	0.8	108	1.2

☞ Severe Wedge
☞ Mild Biconcavity

GE Lunar Report

COMMENTS:



Normal (Grade 0)	<input type="checkbox"/> No Vertebral Deformities Seen		
	Wedge Deformity	Biconcave Deformity	Crush Deformity
Mild (Grade 1)			
Moderate (Grade 2)			
Severe (Grade 3)			

Label	Height (mm)			Percent Deformation		
	Post	Mid	Ant	Wedge	Biconcave	Crush
	Deformity (Grade)					
T4	17.8	18.3	18.2	0.0%	-3.1%	2.5%
T5	19.9	17.9	18.5	6.8%	9.9%	0.0%
T6	19.4	18.1	19.2	1.1%	6.7%	0.0%
T7	21.1	19.6	18.5	12.3%	6.8%	0.0%
T8	20.9	19.4	18.7	10.5%	7.0%	0.0%
T9	21.1	20.2	19.5	7.3%	4.3%	0.0%
T10	22.4	20.2	18.3	18.5%	9.9%	0.0%
T11	24.3	18.8	16.1	33.7%	22.5%	0.0%
T12	26.3	19.4	19.4	26.2%	26.5%	0.0%
L1	28.9	22.5	18.1	37.6%	22.4%	0.0%
L2	27.0	23.7	27.1	0.0%	12.2%	0.5%
L3	27.8	25.4	26.4	5.2%	8.8%	0.0%
L4	27.3	25.1	26.7	2.2%	8.3%	0.0%
Std Dev	1.0	1.0	1.0	5.0%	5.0%	5.0%

Hologic Report

Orthopedic Analysis™

Orthopedic Analysis™

Hip prostheses, metal fastenings, and other artifacts are easily excluded from the analysis region for accurate bone density results.



<http://www.gemedicalsystemseurope.com>



Scan Parameters

Patient Name:	LAW
Patient ID:	
Scan Type:	a SE L.Femur
Scan ID:	A07161309
Scan Length:	15.7 in
Scan Width:	4.5 in
Line Spacing:	0.0095 in
Point Resolution:	0.0355 in
2.40 X 0.01 Coll.	
140/140 kVp	
5.0 mA avg.	
15 seconds	
60 Hz	

**** XRAY DEVICE READY ****
You may begin scanning

Image Not Available

Please position the patient for a SE L.Femur scan

Hologic QDR Workstation

Scan Identification
Patient Name: LAW Scan Type: x Left Hip
Patient ID: Scan ID: A07161307

Positioning SE Femur

Start Position

Scan time: 00:01
Line: 50

For Help, press F1 | Help | 07/16/2013 10:44

Atypical Femur Fracture Scan

SE Femur

This scan is performed when a patient has been on an Osteoporosis therapy for a very long period and the patient is having femur pain.

This is done to detect Atypical Femur Fracture a side effect of bisphosphonates.

The atypical femur fracture is detected by a beaking sign.

If not taken care of a traumatic fracture may occur.



Beaking



OH NO !!!

**Not taken
care of**

SE Femur Image Viewer

The image appears in the SE Femur Image Viewer.

The Viewer allows the operator to view the SE Femur image.

There are controls on the Viewer to change viewing modes, zoom, adjust the brightness and contrast of the image, place Rulers and print a report or image.

Visual Assessment

Visually inspect the image for deformity 1 2 3 particularly along the lateral cortex from the lesser trochanter to the supracondylar flare. Look for focal reaction or thickening along the lateral cortex, which may be accompanied by a transverse radiolucent line. Use the visual tools to increase magnification and adjust the contrast as the changes in the lateral cortex may be subtle.

Select Scan Type:

Whole Body	
AP/Lateral	
Dual-Hip	
AP/Decubitus	
HP Whole Body	
Infant Wholebody	
IVA Imaging	
IVA with BMD	

Patient Name: ADDAMMS, PUGSLEY TH

Patient ID: 1

Use Default Scan Mode

<< Back Next >> Cancel

IVA with BMD




Scan Parameters

Patient Name: Adria, Sophia A
 Patient ID: 9000
 Scan Type: f Lateral
 Scan ID: A1216020G

Scan Length: in
 Scan Width: 5.7 in
 Line Spacing: 0.0397 in
 Point Resolution: 0.0396 in

2.26 X 0.04 Coll.
 140/100 kVp
 2.5 mA avg.
 123 seconds
 60 Hz

**** XRAY DEVICE READY ****
 You may begin scanning

Please position the patient for a Lateral scan




Scan Parameters

Patient Name: Adria, Sophia A
 Patient ID: 9000
 Scan Type: f Lumbar Spine
 Scan ID: A1216020F


Scan Length: in
 Scan Width: 4.5 in
 Line Spacing: 0.0397 in
 Point Resolution: 0.0355 in

2.51 X 0.04 Coll.
 140/100 kVp
 2.5 mA avg.
 41 seconds
 60 Hz

**** XRAY DEVICE READY ****
 You may begin scanning

Please position the patient for a Lumbar Spine scan



The Lateral IVA Scan

The activities involved in performing the Lateral IVA Scan include:

- Confirming Scan parameters for Lateral IVA Scan
- Performing the Lateral IVA Scan
- Using the IVA Image Viewer for the Lateral IVA Scan

Confirming Scan Parameters for Lateral IVA Scan

After the Lateral BMD scan is complete, the scan parameter screen for the Lateral IVA (f SE Lateral Image) will appear. The only scan parameter for the Lateral IVA that can be modified

is the scan length. The scan length is 16.1 inches (40.8 cm).

Hologic QDR Workstation

A10090205 f Lateral

Global ROI
Vertebral Boundaries
Vertebral Bodies
Bone Map
Results

Patient Data
Scan Date: 10/09/2002 15:30
Name: Adria, Sophia A
Pat Id: 9000 Sex: F
Birthdate: 12/30/1947 Age: 54
Height: 64.0 in Weight: 164.0 lb
Ethnic: Hispanic

V Bodies Toolbox
Whole Mode
Line Mode
Point Mode
Select a region
↑ ↓
Include
Exclude
Undo



L2 29 x 25
Dual Energy

Cancel

For Help, press F1

Help 10/09/2002 15:43

Hologic QDR Workstation

A10090205 f Lateral

Global ROI
Vertebral Boundaries
Vertebral Bodies
Bone Map
Results

Bone Map Toolbox

Edit State

- Inactive
- Add Bone
- Delete Bone

Undo
Reset



90 x 106
Dual Energy

Cancel

Patient Data

Scan Date: 10/09/2002 15:30

Name: Adria, Sophia A
Pat Id: 9000 Sex: F
Birthdate: 12/30/1947 Age: 54
Height: 64.0 in Weight: 164.0 lb
Ethnic: Hispanic

For Help, press F1

Help 10/09/2002 15:45

Results button

Hologic Delphi A

A06230303 fLateral

Global ROI
Vertebral Boundaries
Vertebral Bodies
Bone Map
Results

Results Toolbox

BMD
WA-BMD

96 x 106
Dual Energy

Patient Data

Scan Date: 06/23/2003 11:23

Name: Adria, Sophia A
Pat Id: 9000 Sex: F
Birthdate: 12/30/1947 Age: 55
Height: 64.0 in Weight: 164.0 lb
Ethnic: Hispanic

Supine Lateral Analysis

Image not for diagnostic use
k=1.131 d0=37.3 t=9.158 b=1.000
TOTAL BMD CV for L2 - L4 1.0%
C.F. 1.025 1.006 1.000

Region	Area (cm ²)	BMC (g)	BMD (g/cm ²)	T-Score	Z-Score
L2	8.43	7.45	0.884	0.2	2.4
L3	9.93	9.37	0.944	1.1	3.2
L4	11.01	9.48	0.861	0.0	1.5
TOTAL	29.37	26.30	0.896	0.9	2.6

Delphi A SN: 45005
Version 12.0 :3 06/23/2003 11:30

For Help, press F1

Help 06/23/2003 11:30

Results area
Results toolbox

KP1068_002-1233

Orthopedic Analysis™

Orthopedic Analysis™

Hip prostheses, metal fastenings, and other artifacts are easily excluded from the analysis region for accurate bone density results.



<http://www.gemedicalsystemseurope.com>

Trabecular Bone Score (TBS)

TBS iN-sight is a software provided for use as a complement to a DXA analysis. It computes the antero-posterior spine DXA examination file and calculates a score (Trabecular Bone Score - TBS) that is compared to those of the age-matched controls.

TBS is a gray-level textural metric that can be extracted from the two-dimensional lumbar spine dual-energy X-ray absorptiometry (DXA) image.

TBS is related to bone microarchitecture and provides skeletal information that is not captured from the standard bone mineral density (BMD) measurement.

TBS Findings

TBS gives lower values in postmenopausal women and in men with previous fragility fractures than their nonfractured counterparts.




TBS is complementary to data available by lumbar spine DXA measurements.

TBS results are lower in women who have sustained a fragility fracture but in whom DXA does not indicate osteoporosis or even osteopenia.

TBS predicts fracture risk as well as lumbar spine BMD measurements in postmenopausal women.




TBS is associated with fracture risk in individuals with conditions related to reduced bone mass or bone quality.

Healthy Well-structured trabecular bone



TBS = 1.360

Osteoporosis Altered trabecular bone



TBS = 1.100

Risk Class based on minimum hip or spine BMD T-score

		Normal	Osteopenia	Osteoporosis
Risk Class based on Spine TBS	≥ 1.300			
	$1.200 < > 1.300$			
	≤ 1.200			

Adapted from Table 3 in Hans et al. J Bone Miner Res. 2011 Nov;26(11):2762-9


Color coded risk levels for major osteoporotic fracture per 1'000 women per year, based on a $\approx 30'000$ women study.

SPINE TBS REPORT


TBS reference graph

Reference population: USA (NHANES / Medimaps)

TBS L1-L4: 1.259



TBS Mapping




Non diagnostic image

Patient:	Test1, Test1	Patient ID:	
Date of birth:	01/01/1990 - 54,7 years	Acquisition date:	29/09/2014
Height / Weight:	165.1 cm / 60.0 kg	Prescribing doctor:	
Gender / Ethnicity:	Female / White		


SPINE TBS REPORT

TBS reference graph

Reference population: USA



TBS Mapping



Additional results

Region	TBS	BMD	BMD T-Score
L1	1,055	1,207	---
L2	1,257	1,205	---
L3	1,321	1,284	0,0
L4	1,362	1,297	0,0
L1-L4	1,249	1,205	---
L1-L3	1,211	1,252	0,0
L1-L2	1,158	1,235	0,0
L2-L4	1,313	1,252	0,0
L2-L3	1,288	1,273	0,0
L3-L4	1,341	1,291	0,5

Comments

The TBS is derived from the feature of the DXA image and has been shown to be related to bone mineral density and fracture risk. This data provides information independent of BMD values. It is used as a supplement to the data obtained from the DXA analysis and the clinical assessment. The TBS score can assist the health-care professional in assessment of fracture risk and in monitoring the effect of treatments on patients across time. **Overall fracture risk will depend on many additional factors that should be considered before making diagnosis or therapeutic recommendations.** The software does not diagnose disease or recommend treatment regimens. Only the health care professional can make these judgments. USA file: "spinev0701.exe" (TBS analysis done on 25/09/2014, version 3.1.13)
 This DXA system has not been calibrated with a specific TBS phantom. The TBS score has been computed with a generic calibration. These results can be used at the sole discretion of the physician.

TBS can now be applied to FRAX.

A Calculation Tool

Please answer the questions below to calculate the ten year probability of fracture with BMD.

Country: US (Caucasian) Name/ID: About the risk factors

Questionnaire:

1. Age (between 40 and 90 years) or Date of Birth
Age: Y: M: D:

2. Sex Male Female

3. Weight (kg)

4. Height (cm)

5. Previous Fracture No Yes

6. Parent Fractured Hip No Yes

7. Current Smoking No Yes

8. Glucocorticoids No Yes

9. Rheumatoid arthritis No Yes

10. Secondary osteoporosis No Yes

11. Alcohol 3 or more units/day No Yes

12. Femoral neck BMD (g/cm²)
T-Score

BMI: 26.6
The ten year probability of fracture (%)
with BMD

Major osteoporotic	13
Hip Fracture	2.5

If you have a TBS value, click here:

B

FRAX adjusted for TBS

WHO FRAX web site What is TBS? Calculation Tool References TBS web site English

Calculation tool

Country: US (Caucasian)
Name/ID: -
Age: 65
Sex: Female
BMI (kg/m²): 26.6

Please enter the Trabecular Bone Score to compute the ten year probability of fracture adjusted for TBS
Lumbar Spine TBS:

Attention: TBS values are accurate only for patients (women and men) with a BMI in the range [15 – 37 kg/m²]

The 10 year probability of fracture (%)
Adjusted for TBS

Major Osteoporotic Fracture: 17
Hip Fracture: 3.8

00003896
Individuals with fracture risk assessed since April 15, 2015

3D DXA

3D-SHAPER is a software application that registers a 3D statistical model onto the hip DXA scan of the patient and uses a model-based algorithm to create a 3D map of the cortical surface density. It allows clinicians to assess the cortical and trabecular macrostructure in 3D from a standard hip DXA scan.

3D DXA

Incorporates a statistical model built from a database of quantitative computed tomography scans.

Registers the statistical model onto the DXA projection of the patient. The algorithm maximizes the similarity between the projection of the model and the DXA image

Estimates the cortical thickness and density by the fitting a mathematical function onto the density profiles measured along the normal vector at each vertex of the femoral shape

3D-SHAPER
















QCT


3D-SHAPER

QCT

3D-SHAPER

QCT





Grid Plot Viewer


Numeric Data

Output/Global Parameters

13/05/2012 09:11:07

Left Femur Parameter





1. Total (mg)
2. Cortical (mg/cm³)
3. Trabecular (mg/cm³)
4. Cortical Total (mg/cm³)
5. Trabecular Total (mg/cm³)



Exam Date	Age
13/05/2012	69.7

GALGO MEDICAL
 Carreter del Comte d' Urgell 143, Barcelona, Spain


Patient: 2493336 Patient ID: 1001
 Date of birth: 04/05/1942 (70 years) Doctor: Dr. Hippocrates
 Height Weight: 160 cm / 58 kg Acquisition date: 21/02/2012 11:38:40
 Gender: Female Acquisition site: Right Femur
 Ethnic group: White 3D Analysis: 28/05/2012 10:16:41 (02:30)




	Value	T-Score	Z-Score
Cortical sBMD (mg/cm³)	186	-2.7	-1.8
Trabecular vBMD (mg/cm³)	79.2	-3.2	-1.9

General Comments
 [General Comments Here]

Cortical sBMD




Trabecular vBMD



	Value	T-Score	Z-Score
Cortical sBMD (mg/cm ²)	145	-1.0	0.6
Trabecular vBMD (mg/cm ³)	120	-2.0	0.3

Hip Structure Analysis (HSA)



Hip Axis Length (HAL)

Geometric measurements might be used together with densitometric evaluations for a better assessment of hip fracture risk

Measures length and angle of hip axis





Figure 1: Hip image from a Hologic DXA scanner showing positions of thin analysis regions across the femur at the neck (NN region), intertrochanteric (IT) and shaft. On the left are shown typical bone mass profiles used in measurements of geometric properties.



HSA

As HAL increases, fracture risk increases

Each centimeter (10%) increase in Hip Axis length (HAL) increases hip fracture by 50-80% depending on the study




HSA™ Results Summary:

Region	Sub Peri. Width(cm)	Endo Cort. Width(cm)	CSA (cm ²)	CSMI (cm ⁴)	Z (cm ³)	Cort. Thick (cm)	BR
NN	3.84	3.46	3.66	4.31	2.02	0.19	11.2
IT	7.07	5.97	7.89	30.40	7.35	0.55	7.5
FS	3.47	2.12	5.90	6.15	3.44	0.67	2.7
Neck Shaft Angle:	134°						
HAL:	115 mm						


Used forearm protocol.

You can use other protocol to examine other skeletal areas. You just won't get a T-Score



Region	Area (cm ²)	BMC (g)	BMD (g/cm ²)
GLOBAL	23.32	18.50	0.794
R1	11.95	7.65	0.641
Net	11.95	7.65	0.641

Normal Heel



Heel Bone Spur



Sign of Previous trauma to Ankle

