


Scanner Components and Operation – 3




ADVANCED HEALTH EDUCATION CENTER®

EDUCATION ▲ STAFFING ▲ CONSULTING



Major Component -- Gantry

- **Ring-shaped part of the CT scanner**
 - 70-90 cm
 - 15 - 30 degree forward or backward tilt
- **Houses:**
 - Slip rings
 - Generator
 - Cooling system
 - X-ray Tube
 - Filters
 - Collimators
 - Detectors



X-Ray Tubes

- **First and second generation scanners used**
 - fixed anode
 - oil-cooled x-ray tubes
- **Modern CT uses**
 - rotating anode tubes for better tube cooling
 - Focal spot of .5 and 1 mm
 - Small focal spot - better detail
- **Helical scanner tubes rotate for a longer period, thus the tube must be able to sustain high power levels**
- **Many technical and design advances have been made.**





Straton X-ray Tube

- For Multi slice (MSCT)
- Differences from Traditional
 - Entire tube rotates
 - Entire tube immersed in oil
 - Oil is circulated for cooling
 - Deflected Electron beam to 2 focal spots




Due to fastest anode cooling, STRATON never accumulates heat during exposure.





Protective Housing

- X-rays are emitted isotropically, with equal intensity in all directions
- The tube window is the section of tube designed to let x-rays pass through
- These x-rays are called the useful beam
- All others are termed Leakage Radiation



Glass Envelope

- Tube is evacuated and contains two electrodes: cathode and anode
- Most tubes are made of pyrex, some newer tubes have metal construction
- As tungsten filament vaporizes in glass tube arcing can be a problem



Cathode

- Cathode is the negative side of tube
 - Contains: Filament & Focusing Cup.
 - Filament: coil of thoriated tungsten wire and is a source of electrons
 - Tungsten has high melting point 3410°C
 - Focusing cup: filament is embedded in metal focusing cup.
 - Negatively charged focusing cup directs electrons to specific area of anode



Filament

- Dual Filaments: most x-ray tubes have a small and large filament
- Filament provides resistance to flow of electrons and the heat produced causes thermionic emissions to occur.
- These are the electrons that strike the target of the Anode



Anode


Anode is the positive side of x-ray tube

- Two types of anode:
 - stationary (no longer used)
 - Rotating
- Functions of anode
 - Anode target surface is where the high-speed electrons from the filament are suddenly stopped, resulting in the production of X-rays



Focal Spot of Anode

- Small focal spot is associated with small filament.
- Large focal spot is associated with a large filament




Anode Construction

- Anode construction:
 - copper for shaft
 - graphite and molybdenum light weight and dissipate heat.
- Target:
 - area of anode struck by electrons
 - usually constructed of tungsten/rhenium alloy.
- Tungsten's high atomic # (74) results in high efficiency x-ray production.



Target of Anode

- Rotating anode provides approximately 1000 times more area for electron beam to interact than a stationary anode
- Rotating Anode is driven by an electromagnetic induction motor



Rotating Disk/Anode


- Made of an RTM alloy:
 - Rhenium
 - Tungsten
 - Molybdenum
- Small target angle (usually 12 degrees)
- Rotation speed of 10,000 rpm



Line Focus Principle


- Line-focus principle: by angling the target an effective focal spot can be created that is smaller than the actual focal spot.
- When target angle is made smaller, effective focal spot will be smaller.
- X-ray tube target angles vary from between 5° - 15°





X-RAY TUBE

- Line-focus principle simultaneously:
- Improves spatial resolution of the x-ray beam
- As well as the heat capacity of the anode




X-Ray Tube Failure

- Tube failure: tube life is extended by using appropriate radiographic factors of mA, kVp and time.
- Almost all tube failure is related to thermal characteristics of x-ray tube.
- Tube heat must be dissipated through radiation, conduction and convection.
- Tube warm-up procedures must be followed per manufacturers instructions





Gantry Components

What are inside the gantry?



Schematic Representation of the Scanning Geometry of a CT System






Slip-Ring Technology

- **Allows for continuous gantry rotation-faster scan times/no interscan delay**
- **Referred to as continuous rotation, Volumetric or Helical scanning**
- **Fast data collection- required for angiography and continuous acquisition protocols**
- **No high tension cables of the start/stop scanners**



Generator

- Produce high voltage for the creation of x-ray photons
- Increasing the voltage (i.e., kV) will increase the intensity of the x-ray beam. (quality or Hardness of beam)




Generator

- Early scanners used high-voltage stand alone transformer that were bulky with high tension cable from the generator to the x-ray tube
- Now high frequency generator are used
 - Smaller and more efficient at x-ray production
 - Located inside the gantry
- kVp range settings **80, 120 & 140**
 - High kVp reduces patient dose by allowing a lower mAs
 - Reduces heat load by using lower mA setting
- mA range settings **20-500**



Cooling Systems

- The production of x-ray energy results in a large amount of by-product heat
- Blowers, filters, or other devices help to dissipate the excess heat





Filters

- **CT X-ray beams have long and short wavelengths – polychromatic.**
 - **Original scanner had monochromatic beam that satisfied the Lambert-Beer exponential attenuation law**
 - **Beam needs to be monochromatic or nearly monochromatic for reconstruction purposes. (Lambert-Beer exponential attenuation law)**
- **CT Filter has dual Purpose**
 - **Removes long wavelengths (low energy) x-rays that contribute to patient dose but DO NOT help form the image**
 - **This is called beam hardening**
 - **Filtration shapes the energy distribution across the radiation beam to produce uniform beam hardening. This reduces beam hardening artifacts**




Bowtie Filters and Beam Hardening



Collimation

- Restricts the x-ray beam to a specific area
- Reduces scatter radiation
- Improves contrast resolution
- Decreases patient dose



Collimation

- **Two collimators in CT:**
 - **Pre-patient**
 - **Post patient or Pre-detector**
- **Pre-patient collimator**
 - **positioned near the x-ray source**
 - **Controls the patient dose**
- **Post Patient Collimator**
 - **Removes scatter and improves image quality**



Collimators

- **Both sets of collimators insure a constant beam width**
- **Improves Contrast Resolution**
 - Contrast Resolution is the ability to differentiate small density differences on the image
 - Reducing scatter always improves contrast resolution
- **Post patient collimators control the slice thickness by widening or narrowing the x-ray beam.**
- **As the slice thickness changes so does voxel size.**



Detector - Scintillation (solid state)

- **Converts x-ray energy into light, and then the light is converted into electrical energy**
- **Consist of scintillation crystals coupled to a photomultiplier (PM) tube**
- **Early scanners used sodium iodide crystals couple to PM tubes**
- **Today most scintillators are cadmium tungstate or high-purity rare earth oxides**
- **The detector efficiency rate of solid state detectors is 99%.**
- **Automatically recalibrated between scans.**



Detector Gas - Ionization Detectors

- **Based on the principle of ionization**
- **Gas-filled detector array uses small ion chambers filled with high-pressure xenon or other gas**
- **Each ion chamber is about 1mm wide with essentially no inter-space**
- **Has excellent stability and fast response times but detector efficiency is only 50%.**



Detector Characteristics

- **Efficiency :**
 - The ability to capture, absorb, and convert x-ray photons to electrical signals.
 - **Capture efficiency:**
 - The efficiency with which the detectors can obtain photons transmitted thru the patient
 - Size of detector face
 - Spacing of detectors
 - **Absorption efficiency:**
 - the number of photons absorbed by the detector
- **Stability**
 - Steadiness or ability to stay in calibration
- **Response Time**
 - Speed of detecting and recovering to capture a new photon



Data Acquisition System (DAS)

- **Positioned Between**
 - **Detector array and Computer**
- **3 functions:**
 - **Measures the photons that pass through the patient and strike the detectors**
 - **Converts data from analog to digital**
 - **Transmits binary data to computer**

Data Acquisition System

- The DAS consists of the following parts




- X-ray photons come on the detector.
- The detector detects the intensity in form of current.
- The current is converted into voltage.
- The analog integrator removes spikes.
- The analog signal is converted into digital form.
- This signal can now be processed and reconstructed in the computer.



Summary of Process

- Tube
- Filter
- Collimators
- Patient
- Post patient collimators (not shown)
- Detector
- DAS
- Computer to make image



Patient Table


- The process of moving the table within the gantry by a specified measure is called incrementation, feed, step, or index
- The degree to which a table can move horizontally is called the scanable range
 - This determines the extent a patient can be scanned without repositioning



Consoles

CT Consoles - 2 types

- **Operating**
 - Key Boards
 - Mouse
 - Flats Screens for viewing and programing
- **Viewing**




Archival Devices

- Old Systems
 - Magnetic Tapes
 - Magnetic
 - Floppy and Hard
 - Film
- New Systems
 - CDs
 - PACS (Picture Archiving and Communicating System)



PACS



ADVANCED HEALTH EDUCATION CENTER®

EDUCATION ▲ STAFFING ▲ CONSULTING

